VitiCELL®

Kit for preparation of a cell suspension

Kit for preparation of a cell suspension

INSTRUCTIONS FOR USE

DESCRIPTION

Viticell® is a single-use, stand-alone device for performing a cell suspension. Viticell® is a kit comprised of reagents and attached instruments: enzymatic solutions, sterile application solutions and devices. The kit also includes a portable instant heat pack containing an aqueous solution saturated with sodium acetate, which provides heat through crystallisation.

Starting with a thin-skin biopsy, this device enables, through disintegration of the graft, a cell suspension made up of a mixed population, primarily of keratinocyte basal cells but also containing Langerhans cells, melanocytes and fibroblasts.

The suspension obtained allows culturing of the surface of a wound on the same patient.

The Viticell® kit must be utilised according to a physician's recommendations. The sampling, suspension and therapeutic application of autologous cells must be performed by a physician in observance of current laws.

INDICATIONS AND DOSAGE

Viticell® is a device that enables cell suspension and, in particular, it can be dedicated to skin repigmentation, for example as part of non-progressive vitiligo or hypochromatic post-traumatic scars.

Viticell® contributes the reagents and devices required for preparing an autologous skin cell suspension from a biopsy measuring 4 to 10 cm² of thin skin (0.2-0.3 mm). The material supplied also allows drop-by-drop application onto a wound with a surface area measuring up to 10 times the size of the biopsy. The surface area to be treated must not exceed this ratio of 10; the recommended ratio is 5 times the size of the biopsy.

CONTRAINDICATIONS

The cell suspension must not be applied to an infected or necrotic wound.

Viticell® must not be utilised with patients who are hypersensitive to hyaluronic acid or rtrypsin.

PRECAUTIONS FOR USE

Viticell® allows a cell suspension to be obtained for exclusively autologous use (the treatment must be applied to the patient from whom the biopsy was taken).

Viticell® is for a single use; its components must be utilised only for Viticell® and are not reusable in order to avoid any risk of infection to the patient.

All solutions supplied in the Viticell® kit that are intended to come into contact with the cells are sterile.

For handling the kit, follow the recommended precautions for use.

The kit's components must not be frozen or re-sterilised. The kit's components must not be mixed with other products.

The enzyme utilised for tissue separation (rtrypsin) is recombinant in origin. The safety controls implemented by the enzyme's manufacturer enable minimisation of pathogenic contamination risks even if these risks cannot be totally ruled out.

Handling of the Viticell® kit takes place in a sterile field. Once opened, Viticell® must be utilised immediately and discarded after use.

Before beginning to use Viticell®, check that all of the Viticell® components are present.

Check the completeness of the heat pack in the box. The heat pack must not have been activated (the heat pack appears hard and opaque white in colour when it has been activated).

In order to guarantee the optimal viability of the cells, it is recommended to process the biopsy immediately after sampling and to perform the graft on the lesion to be treated immediately after this process

Do not utilise antiseptics containing quaternary ammonium salts on the skin to be treated (i.e. Biseptine®, Mercryl®).

Keep out of reach of children.

DO NOT INJECT.

The benefit of using the kit for repigmentation in pregnant and breast-feeding women is to be assessed by the Physician.

To store the kit, see the STORAGE section.

ADVERSE EFFECTS

Possible side effects can be rarely caused by the cell transplantation technique at the donor site or on the treated zone: local infection, pigmentation disorders, scars, dryness, pruritus, pain, erythema.

Any adverse effect associated with Viticell® must be reported immediately to the manufacturer.

(contact: vigilances.fr@ibsa-group.com).

KIT COMPONENTS

Viticell® consists of:	
1 x 10-ml pre-filled syringe containing 6 ml of digestion medium (PBS) (DB)	1. <i>Digesting Buffer</i>
1 x 10-ml pre-filled syringe containing 5 ml of rinsing medium (PBS) (WB)	2. <i>Washing Buffer</i>
1 x 2.25-ml pre-filled syringe containing 1.5 ml of preparation medium (PBS) (CB)	3. <i>Collecting Buffer</i>
1 x 2.25-ml pre-filled syringe containing 1 ml of hyaluronic acid 2% for cell suspension and application (HA)	4. <i>Hyaluronic Acid 2 %</i>
1 x 5-ml amber bottle containing 7 mg of rtrypsin (DC)	<i>Digesting Cocktail</i>
1 x instant heat pack	
1 x compartmentalised Petri dish	
1 x empty container	
3 x 18G transfer needles	

Sterilisation of components:

The syringes containing the sterile digestion medium (DB 1. *Digesting Buffer*), the sterile rinsing medium (2. *Washing Buffer*) and the sterile preparation medium (3. *Collecting Buffer*) were sterilised using ethylene oxide, and then filled and steam-sterilised.

The 2% hyaluronic acid solution (HA 4. *Hyaluronic Acid 2 %*) was steam-sterilised.

The rtrypsin (DC *Digesting Cocktail*) was sterilised by irradiation.

The compartmentalised Petri dish was sterilised by irradiation.
The empty upright container was sterilised by irradiation.
The 18G transfer needles were sterilised by ethylene oxide.

STORAGE

Upon receipt of the Viticell®, check that the packaging is intact.
The Viticell® kit must be stored between +15°C and +25°C. The kit's components must be stored in their packaging until use; the enzyme (DC *Digesting Cocktail*) is particularly sensitive to light, so it must remain in its original bottle. Do not use Viticell® or its components if the use-by date indicated on the packaging is expired. Do not use Viticell® or its components if the packaging is damaged or if, upon visual inspection, it has been opened (check tamper-proof safety label); if necessary, contact IBSA Pharma SAS immediately.

DISPOSAL

Viticell®'s components must be discarded after a single use in appropriate biomedical waste receptacles.

REQUIRED MATERIALS (NOT SUPPLIED)

Appropriate sterile surgical field;
Sterile gloves;
Sterile forceps;
Sterile scalpel;
Protective goggles and clothing;
Surgical mask;
Instruments required for biopsy, scalpel or dermatome as well as associated materials (sterile knife and forceps), antiseptic solution and sterile saline solution;
Instruments required for dermabrasion (laser or equivalent technique) and associated materials;
Primary and secondary bandages.

SUMMARY OF THE ADMINISTRATION METHOD

The process takes about 1 hour (see details in the section titled INSTRUCTIONS FOR USING VITICELL®):

1. The reagents are prepared (5 minutes).
2. A thin-skin sample measuring 4 to 10 cm² is taken on a healthy skin area (5 min).
3. A cell digestion / isolation phase is performed thanks to the device, in 2 stages (20 min):
 - Tissue digestion (15 min),
 - Dermoepidermal separation (5 min).
4. A suspension phase of the cells isolated in hyaluronic acid (5 min).
5. A culturing phase on the prepared lesion (5 min).

Do not inject the cell suspension.

A primary bandage must be placed on the area to maintain the graft.

In the case of partial depigmentation, a second graft may be considered in the months following the first application.

INSTRUCTIONS FOR USING VITICELL®

It is recommended that you read through all instructions once before starting to use the kit.

Prepare the patient's sampling area beforehand: if necessary, disinfect it and use local anaesthetic.

1. PREPARING THE DEVICE

Prepare a sterile surgical field.


Take the three trays out of the box and place them on the sterile field. The heat pack must stay on its tray; this will allow positioning of the Petri dish and optimisation of the heating. Unwrap the Petri dish and place it on the heat pack. The empty upright container must be positioned in the location intended for this purpose to the right of the heat pack.

THE NEEDLES ARE FOR TRANSFERRING FLUID AND PREPARING THE SUSPENSION; THEY ARE NOT INTENDED FOR INJECTING THE PRODUCT INTO THE PATIENT. DO NOT INJECT EITHER THE CELL SUSPENSION OR ANY OTHER KIT COMPONENT INTO THE PATIENT.

2. PREPARING AND HEATING THE ENZYME SOLUTION

- Open the boxes containing the syringes.
- Remove the bottle's membrane seal and stopper (DC *Digesting Cocktail*).
- Wearing sterile gloves and a mask, screw onto the syringe (DB 1. *Digesting Buffer*) a transfer needle.
- Inject the contents of the pre-filled syringe (DB 1. *Digesting Buffer*) into the bottle of powder (DC *Digesting Cocktail*).
- Homogenise the solution and, using the syringe, transfer it into one of the two compartments of the Petri dish.

3. TAKING A SKIN SAMPLE

Activate the instant heat pack by pressing on the metal plate; the crystal solidifies instantly, the pack begins to heat up and it turns an opalescent white colour.

Put the Petri dish on top of the activated heat pack and wait 5 min, during which time you can perform the biopsy.

N.B.: Be sure not to exceed 5 minutes, since the heat pack has a limited heating duration.

A thin-skin biopsy measuring 0.2 – 0.3 mm thick and 4 to 10 cm² in surface area must be taken from a donor site with pigmentation and texture that are similar to those of the site to be treated. Determine the biopsy's size according to the surface area to be treated. The surface area to be treated must not exceed this ratio of 10; the recommended ratio is 5 times the size of the biopsy. Clean the donor site with an antiseptic solution, then rinse it with a sterile saline solution. Local anaesthetic can be utilised, if necessary, subcutaneously.

Use of a dermatome is recommended to obtain a sufficiently thin biopsy (0.2 – 0.3 mm thick) in order to facilitate its digestion and promote healing of the donor site without scarring. Next, treat the donor site as usual.

4. ENZYMATIC DIGESTION OF THE BIOPSY

With the help of sterile forceps (not supplied), transfer the biopsy into the compartment of the Petri dish that contains the *Digesting cocktail* solution that is still on the heat pack, and close it using the cover. Heat the Petri dish containing the biopsy for **15 minutes** on the instant heat pack.

5. RINSING THE BIOPSY

After the 15 minutes, deposit the contents of the pre-filled syringe (**WB 2. Washing Buffer**) into the second compartment of the Petri dish (without installing a needle).

With sterile forceps (not supplied), transfer the biopsy into this second compartment for a quick rinsing in the (**WB 2. Washing Buffer**), then remove the biopsy and place it inside the Petri dish's cover while being careful to keep the dermoepidermal junction pointing upward.

Using sterile forceps and scalpel (not supplied in the kit), separate the dermal part from the epidermal part.


6. COLLECTING SKIN CELLS

Deposit the preparation medium from the pre-filled syringe (**CB 3. Collecting Buffer**) onto the biopsy fragments, drop by drop, without installing a transfer needle. With the scalpel, repeatedly scrape the cells from the junction surface and cut the epidermal part into small fragments to achieve a mixture of homogeneous cells.

Open the empty upright container and position it in the location intended for this purpose to the right of the heat pack.

Screw the transfer needle onto the empty syringe (CB 3. *Collecting Buffer*) and take out all of the cell suspension by tilting the container slightly.

Replace the stopper on the bottle and leave at room temperature during the preparation time of the site to be treated.

7. PREPARING THE SITE TO BE TREATED

The site to be treated must undergo dermabrasion using a laser or an equivalent technique. The area must be disinfected and anaesthetised beforehand, if necessary.

Do not utilise antiseptics containing quaternary ammonium salts on the skin to be treated (i.e. Biseptine®, Mercryl®).

8. PERFORMING THE CELL SUSPENSION

Install the transfer needle onto the pre-filled syringe containing the hyaluronic acid solution (HA 4. *Hyaluronic acid 2.0%*) and inject this into the container that holds the cell solution.

Homogenise the cell suspension by going back and forth with the syringe

(HA 4. *Hyaluronic acid*).

Next, remove the suspension in the syringe twice (the volume of the syringe does not enable the full suspension to be removed).

Apply the cell suspension immediately.

9. APPLYING THE CELL SUSPENSION

Using the syringe (HA 4. *Hyaluronic acid*) with the transfer needle still installed, deposit the suspension drop by drop onto the area to be treated.

Do not inject the cell suspension.

10. POSTOPERATIVE TREATMENT

Apply a bandage to maintain a moist environment. The bandage must have fine pores, be non-adhesive and have low absorption; it needs to stay on for a week. A secondary bandage can be applied.

The follow-up must comply with the standard protocols for skin grafts.

Phototherapy may be considered after a cell suspension graft using Viticell® in accordance with the usual protocols for treating vitiligo.

PROBLEMS YOU MAY HAVE**The heat pack was already activated.**

You can reactivate the heat pack by immersing it in boiling water for 5 to 10 minutes, then letting it cool down.

The heat pack broke inside the box.

Do not start the process. Instead, contact the supplier and ask them to provide you with a new heat pack.

The digestion cocktail powder (rtrypsin) does not dissolve completely.

Be sure to resuspend the powder (rtrypsin) completely with the digestion medium by using the syringe provided or by turning the contents of the digestion cocktail bottle upside down several times. It is possible that a small part of the rtrypsin does not dissolve in the reconstituted solution, but the process can continue.

The dermoepidermal separation is difficult to perform.

Check that the heat pack was activated properly for at least 15 minutes during the enzymatic digestion phase. If the biopsy is too large or too thick, the dermoepidermal separation will take longer. Be sure to observe the size and thickness criteria of the biopsy to be taken.

Missing components in the Viticell® box.

Despite all the verification controls implemented within the manufacturing process, if a component is missing in the Viticell® box, contact the manufacturer.

Incident: the biopsy comes in contact with an unsterile area.

There is a risk that the biopsy will be contaminated; do not continue preparing the cell suspension.

Aspiration problems with the suspension due to large aggregates.

The presence of aggregates in the cell suspension is normal and possible, with no consequence in the progress of the process. However, if larger-sized aggregates were to impede the flow inside the syringe, reject those and take samples only from the homogenous cell suspension.

IF YOU HAVE PROBLEMS, DO NOT HESITATE TO CONTACT THE MANUFACTURER.

Manufacturer: IBSA Pharma SAS
Parc de Sophia-Antipolis - Les 3 Moulins
280, rue de Goa - 06600 ANTIBES
France
Tel: +33 (0)4 92 91 15 60.
Fax: +33 (0)4 92 91 15 30
www.ibsa-pharma.fr
Email : mail.fr@ibsa-group.com

Date of revision of the leaflet: 2021-10

	18G Transfer needles EC Marking BDTM Blunt Fill Needles Manufacturer: Becton Dickinson and Company 1 Becton Drive - Franklin Lakes - NJ 07417 USA	 STERILE EO
--	---	-------------------------------

Meaning of the symbols

Consult instructions for use.

Sterilized using ethylene oxide (needles and empty syringes)

Caution: read the warnings carefully

Sterilized using irradiation (DC + petri dish + empty container)

Do not re-use

Sterilized using steam (HA solution + DB + WB + CB)

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Fabbricante: IBSA Pharma SAS