
Aliaxin[®] EV

ESSENTIAL VOLUME

Cross-Linked
Hyaluronic Acid

Hydro**lift**[®]
ACTION

hydro**ACTION** ●○○○○○
lift**ACTION** ●●●●●●
CROSS-LINKING ■■■■■■


1


2


3

4

5

IT

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EN

ASSEMBLING THE NEEDLE ON THE SYRINGE

For optimal handling of Aliaxin® EV Essential Volume products, it is important to assemble the needle on the syringe in accordance with the 5 steps shown in the diagram (see figures 1 to 5).

1

Firmly hold the glass syringe and the luer-lock between the thumb and the index fingers. Grasp the protective cap and unscrew it.

2

Remove needle cap, then firmly insert the needle thread on the end of the syringe.

3

Gently screw the needle in a clockwise direction. Continue to screw until the thread comes into contact with the body of the syringe.

4

Check that the needle is correctly fitted onto the thread. Failure to follow these precautions could lead to risks of the needle dropping out and/or leaks from the luer-lock.

5

With the syringe held in one hand and the cap in the other, remove the needle protective device by pulling firmly.

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The product is only to be used by medical staff as described in this information leaflet. Do not use for applications different than those indicated in this leaflet.

DESCRIPTION

Aliaxin® EV Essential Volume is a reabsorbable medical device (physiological, apyrogenous sterile gel) to be used as a filler for the correction of deep facial cutaneous sagging and to increase volume. The principle component is cross-linked Hyaluronic Acid of non-animal origin, produced by bacterial fermentation.

COMPOSITION

Cross-linked hyaluronic acid.....25 mg/g

Phosphate buffered saline, Water for injectable solutions.....1 g

PACKAGE

Aliaxin® EV Essential Volume is available in a 2ml package. Each package contains:

- package insert leaflet
- 2 sealed blisters, each containing 1 sterile pre-filled mono-use syringe of 1 ml
- Adhesive labels to stick to the patient's medical record for product traceability
- 4 sterile needles (27G X 13mm TSK C € 0123 - 27G X 19mm TSK C € 0123)

MODALITY OF ACTION

Aliaxin® EV Essential Volume must be injected into the sub-cutaneous tissue for the supplementation of the extra-cellular matrix and the volumetric increase of the tissue, that is, for the treatment of cutaneous depressions caused by wrinkles, scarring and volume loss.

AREAS OF USE

Aliaxin® EV Essential Volume is an injectable medical device conforming to the Dir. 93/42/CEE MOD, indicated for the treatment of imperfections in the following areas of the face:

- nasolabial folds
 - post-traumatic and acne scars
 - areas of the face which require enrichment of the facial tissue (cheeks, chin, cheekbones) with temporary volumetric increase through growth of the soft tissues
- The outcome of the treatment depends on the type of skin and the nature of the imperfection: the less evident the imperfection, the better the result.

INSTRUCTIONS FOR USE

Before any form of intervention with Aliaxin® EV Essential Volume, the doctor must carry out an adequate case history and a comprehensive assessment of the conditions of the patient, in order to guarantee the total absence of contraindications to the implant.

The areas to be treated must be identified and evaluated bearing in mind the criteria relative to lines and symmetries to follow.

Local anaesthetic may be used in order to guarantee the necessary comfort to the patient.

The doctor must inform the patient in advance of the modalities of the intervention, its nature, warnings, precautions and possible individual outcomes, potential adverse responses, the expected duration of the treatment and the eventuality of a supplementary intervention for the maintenance and/or definition in detail of the achieved result.

The area of intervention must be cleansed with antiseptic solution before proceeding to the implant.

Extract the syringe from the blister pack, unscrew the top and assemble the needle together with the protection. The latter is only to be removed right before the intervention. Aliaxin® EV Essential Volume is administered with a sterile needle (27G X 13mm - 27G X 19mm) in compliance with the Luer-Lock standards with normalised connections, which is included in the pack.

The injection must be carried out in correspondence to the sub-cutaneous superficial layer; the procedure is, however, at the discretion of the doctor and depends on the correction desired and the method of implantation adopted.

At the end of the procedure it is advisable to carry out a delicate massage of the treated area in order to allow for better distribution of the product.

DOSAGE AND ADMINISTRATION

The volumetric grading printed on the syringe has an indicative value: it is left to the doctor's discretion to choose the dosage to be used for each individual case. The intervals at which the treatment should be repeated depend on various factors, regarding both the physiology of the patient (type of skin, individual metabolism, anatomy, age) and their lifestyle. Another element to consider is relative to the implant techniques used. For maintenance of the obtained results it is advisable to repeat the treatment periodically, every 8-12 months.

PRECAUTIONS

During the treatment the usual precautions in the case of percutaneous procedures must be adopted. The risks are the common ones of infection relative to the type of treatment.

Aliaxin® EV Essential Volume must not be used on patients affected by:

- infective or inflammatory processes near the area of intervention
- ongoing cutaneous allergies
- known hypersensitivity to cheloids
- allergy to the components
- immune system illnesses
- serious pathological skin conditions
- problems with coagulation or anti-coagulating therapies in progress

Close to treatment the patient should avoid consumption of substances which act on blood fluidity (aspirin, NSAID, Vitamin E), in order to minimise the possibility of bruising or bleeding of the injected areas.

Aliaxin® EV Essential Volume may not be implemented in areas where there are mammary, tendon, bone or muscular implants present.

Following the treatment, and until the complete absorption of swelling and reddening, the implanted areas must not be exposed to make-up, excessive heat (sun, UV tanning sessions or laser) or intense cold.

After use, needles and syringes must be disposed of according to the procedures indicated for hospital waste.

ADVERSE REACTIONS, WARNINGS AND CONTRAINDICATIONS

As with any percutaneous injection, there can be incidents of:

- inflammatory reactions (reddening, edema, etc), sometimes associated with itching and pain to the touch
- hematoma
- hardening or nodules at the point of injection
- colouring or discolouring of the skin at the point of injection

These reactions usually disappear within a few days. If they persist for more than a week the patient should contact their doctor immediately.

Aliaxin® EV Essential Volume is an intra cutaneous injection and must not be injected into blood vessels.

There are no reports of problems due to excessive doses or interaction with other medicines.

Do not use in pregnancy.

Aliaxin® EV Essential Volume is available in single dose sterile packages.

The product must not be reused for future applications or on other patients.

Once open the product must be used immediately.

Any remaining product must be disposed of.

Do not use the product if the packaging is damaged.

Do not mix with other injectable substances and do not carry out other implants together with Aliaxin® EV Essential Volume.

The implant must be carried out in an aseptic environment observing the appropriate techniques.

Fill in the adhesive label contained in the box and stick it to the patient's medical records kept at the doctor's office.

Keep out of the reach of children.

STORAGE

Aliaxin® EV Essential Volume must be stored at a temperature between 2° and 28°C. Do not freeze. Do not expose to heat.

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Emergo Europe B.V.
Rose Pharma S.A.
IBSA Farmaceutici Italia Srl